

Environmental DNA Metabarcoding Sequencing Package

—Accurate and powerful, providing an ideal tool for environmental monitoring and scientific research

Features

- One-stop product solution**
 Integrated sample treatment, library preparation, sequencing, and analysis, a true one-stop product solution.
- Wide range of species covered**
 The monitored species cover bacteria, fungi, algae, plants, planktonic animals, benthic animals, and fish, meeting various application needs.
- Excellent data quality**
 The unique DNBSEQ technology provides high-quality sequencing data.
- Automation friendly**
 Combined with automated nucleic acid extraction instruments, automated sample preparation systems and analysis software, to achieve automated sample preparation and data analysis, save labor.
- Powerful analysis function**
 Self-developed software provides various analysis functions, such as OTU analysis, species composition analysis, alpha diversity analysis, beta diversity analysis, biomarker analysis, and more.
- Support large-scale parallel sequencing**
 Up to 4608 barcodes can be provided for library preparation, meeting the demands of ultra-high-throughput sequencing.

Introduction

Strengthening environmental monitoring is crucial for environmental protection. Environmental biomonitoring utilizes the responses of biological individuals, populations, or communities to assess the environmental conditions. It can be applied to monitor atmospheric, aquatic, and soil environments, serving as an important approach to evaluate natural environmental conditions and biodiversity, and is also a significant aspect of environmental research. Environmental DNA sequencing technology is one of the most revolutionary technologies in the field of ecological and environmental sciences since the 20th century. It involves extracting DNA from environmental media such as water, soil, and sediments, followed by PCR amplification and high-throughput sequencing of specific DNA fragments in the genome, enabling qualitative and quantitative analysis of the biological community.

MGI Environmental DNA metabarcoding sequencing package is a combination of self-developed reagents for nucleic acid extraction, ATOplex multiplex PCR library preparation, MGISP-960 automated sample preparation system, DNBSEQ-G99 and DNBSEQ-E25 sequencers, and MetaSIS analysis software. This comprehensive solution covers the entire process from sample processing to analysis reporting. It enables rapid and accurate analysis of the biological community structure in environmental samples, providing tool for ecological environmental monitoring and environmental scientific research.

Table 1. Product parameters

| Panel | Applicable species types | Recommended read length | Recommended data per sample | Sequencer | Recommended sample number /FC |
|----------|------------------------------|-------------------------|-----------------------------|----------------|-------------------------------|
| MiFish | fish | PE150 | ≥65K reads | DNBSEQ-G99/E25 | 500/72 |
| Ac12S | fish | PE300 | | | |
| 16S V3V4 | bacteria | | | | |
| COI | benthic animals, zooplankton | PE300 | DNBSEQ-G99 | 500 | |
| 18S V4 | algae, invertebrates | | | | |
| ITS1 | Fungus | | | | |




Figure1 Workflow

Performance

The following three sets of data were to evaluate the performance of the solution, such as repeatability, species detection rate, and quantitative assessment of species abundance, using simulated samples with known species composition and abundance.

Samples S1-S3 are three replicates of mixed samples containing 8 bacterial DNA standards. The samples were prepared to 16S V3V4 library, sequenced and analyzed on the DNBSEQ-G99. The results are shown in Table 2 and Figure 2. All expected species were detected, and the measured abundance similar with the theoretical abundance. The three replicates showed minimal deviation.

Table 2. Species abundance analysis results

| Species | Theoretical abundance | S1 Measured abundance | S2 Measured abundance | S3 Measured abundance |
|-----------------------|-----------------------|-----------------------|-----------------------|-----------------------|
| <i>Pseudomonas</i> | 4.20% | 6.83% | 7.86% | 6.94% |
| <i>Escherichia</i> | 10.10% | 12.62% | 9.15% | 13.29% |
| <i>Salmonella</i> | 10.40% | 11.48% | 8.12% | 9.97% |
| <i>Lactobacillus</i> | 18.40% | 16.24% | 20.45% | 16.58% |
| <i>Enterococcus</i> | 9.90% | 10.63% | 11.14% | 13.02% |
| <i>Staphylococcus</i> | 15.50% | 12.44% | 18.34% | 12.58% |
| <i>Listeria</i> | 14.10% | 15.40% | 12.56% | 13.97% |
| <i>Bacillus</i> | 17.40% | 14.36% | 12.37% | 13.64% |




Figure 2. Species composition and abundance stacked histogram

Samples S4-S6 are three replicates of mixed samples containing 10 species fish tissue. The samples were prepared to MiFish library, sequenced on the DNBSEQ-E25 platform and analysis. The results are shown in Table 3 and Figure 3. From the results, we can see all expected species were detected, and the measured abundance similar with the theoretical abundance. Additionally, the three replicates showed excellent reproducibility.

Table 3. Species abundance analysis results

| Species | theoretical abundance | S4 measured abundance | S5 measured abundance | S6 measured abundance |
|---|-----------------------|-----------------------|-----------------------|-----------------------|
| <i>Acanthopagrus</i> | 15% | 21% | 21% | 21% |
| <i>Carassius cuvieri</i> | 3% | 2% | 2% | 2% |
| <i>Channa maculata</i> | 6% | 8% | 8% | 8% |
| <i>Ctenopharyngodon idella</i> | 10% | 12% | 12% | 13% |
| <i>Hypophthalmichthys nobilis /Hypophthalmichthys molitrix*</i> | 5% | 4% | 4% | 4% |
| <i>Larimichthys crocea</i> | 16% | 14% | 15% | 15% |
| <i>Lateolabrax japonicus /Laterolabrax maculatus**</i> | 17% | 15% | 15% | 14% |
| <i>Oreochromis niloticus</i> | 5% | 3% | 3% | 3% |
| <i>Sciaenops ocellatus</i> | 12% | 9% | 9% | 9% |
| <i>Sebastiscus marmoratus</i> | 12% | 12% | 12% | 12% |




Figure 3. Species composition and abundance stacked histogram

Samples S7 and S8 are mixed samples containing 10 fish tissue species but with different proportions. The samples were prepared to AC12S library, sequenced and analyzed on the DNBSEQ-G99. The results are shown in Table 4 and Figure 4. We can see all expected species were detected, and the measured abundance similar with the theoretical abundance.

Table 4. Species abundance analysis results

| Species | S7 theoretical abundance | S7 measured abundance | S8 theoretical abundance | S8 measured abundance |
|---|--------------------------|-----------------------|--------------------------|-----------------------|
| <i>Acanthopagrus</i> | 3.79% | 6.81% | 1.94% | 1.76% |
| <i>Carassius cuvieri</i> | 21.00% | 21.97% | 17.19% | 19.01% |
| <i>Channa maculata</i> | 14.73% | 19.70% | 15.92% | 23.68% |
| <i>Ctenopharyngodon idella</i> | 10.32% | 8.65% | 13.03% | 13.35% |
| <i>Hypophthalmichthys nobilis /Hypophthalmichthys molitrix*</i> | 11.47% | 14.86% | 9.88% | 14.36% |
| <i>Larimichthys crocea</i> | 4.01% | 4.30% | 2.01% | 1.31% |
| <i>Lateolabrax maculatus</i> | 4.24% | 3.48% | 6.13% | 5.41% |
| <i>Oreochromis niloticus</i> | 18.49% | 11.52% | 15.98% | 8.49% |
| <i>Sciaenops ocellatus</i> | 3.97% | 2.32% | 5.95% | 3.39% |
| <i>Sebastiscus marmoratus</i> | 7.99% | 6.39% | 11.97% | 9.25% |




Figure 4. Species composition and abundance stacked histogram

Real environmental samples mainly include water, sediment, and soil. Water are suitable for monitoring fish, zooplankton, phytoplankton, and benthic organisms. Sediment are suitable for monitoring benthic organisms and fish. Soil contain a large number of microbial communities. In the following tests, we detected fish in water, benthic organisms in sediment, and fungi in soil.

Samples D3-D8 are water samples collected from 6 different locations in Donghu Lake, Wuhan. Each location was sampled three times. After filtration, the filter membrane were processed using MGI's solution for nucleic acid extraction, MiFish library preparation, DNBSEQ-G99 sequencing and analysis. The results of quality control and data processing are shown in the table 5. We can see the Q30 > 94% for all samples, indicating excellent sequencing quality. Each sample was trimmed to 50,000 reads and conducted to low-quality filtering, denoising, merging, and chimeric removal to generate feature sequences, also known as OTUs.

Table 5. Quality control and data processing results

| Sample ID | Raw reads | Target reads | Down sampled reads | Filtered reads | Q30(%) | Denosed reads | Merged reads | Chimeric reads | Feature reads |
|-----------|-----------|--------------|--------------------|----------------|--------|---------------|--------------|----------------|---------------|
| D3-1 | 815929 | 58832 | 50000 | 47434 | 96.17 | 47220 | 43499 | 1961 | 41538 |
| D3-2 | 831507 | 58949 | 50000 | 47327 | 96.09 | 47104 | 42609 | 1860 | 40749 |
| D3-3 | 845240 | 59063 | 50000 | 47385 | 95.71 | 47106 | 44854 | 3223 | 41631 |
| D4-1 | 856348 | 58529 | 50000 | 47643 | 96.77 | 47333 | 41628 | 1471 | 40157 |
| D4-2 | 734924 | 58997 | 50000 | 47638 | 96.71 | 47486 | 45974 | 2068 | 43906 |
| D4-3 | 856431 | 58840 | 50000 | 47626 | 96.99 | 47352 | 44918 | 2304 | 42614 |
| D5-1 | 749180 | 59064 | 50000 | 47522 | 96.03 | 47258 | 42392 | 2887 | 39505 |
| D5-2 | 724429 | 58958 | 50000 | 47553 | 96.46 | 47352 | 42476 | 2329 | 40147 |
| D5-3 | 740403 | 58817 | 50000 | 47454 | 96.34 | 47225 | 41429 | 2633 | 38796 |
| D6-1 | 711722 | 58710 | 50000 | 46795 | 94.86 | 46558 | 42725 | 4450 | 38275 |
| D6-2 | 683555 | 58926 | 50000 | 47461 | 96.3 | 47213 | 43249 | 5111 | 38138 |
| D6-3 | 667056 | 58724 | 50000 | 47513 | 96.43 | 47189 | 42167 | 5503 | 36664 |
| D7-1 | 1136152 | 58868 | 50000 | 47186 | 95.62 | 46946 | 43338 | 1567 | 41771 |
| D7-2 | 673020 | 59033 | 50000 | 47415 | 96.62 | 47231 | 44439 | 2473 | 41966 |
| D7-3 | 989925 | 58825 | 50000 | 47087 | 95.2 | 46687 | 44323 | 2607 | 41716 |
| D8-1 | 895021 | 58825 | 50000 | 47236 | 95.67 | 47036 | 45033 | 4283 | 40750 |
| D8-2 | 862580 | 58705 | 50000 | 47422 | 95.94 | 47122 | 44556 | 2717 | 41839 |
| D8-3 | 631941 | 58872 | 50000 | 47487 | 96.13 | 47208 | 44753 | 2410 | 42343 |

After obtaining the OTUs, the system aligns the OTU sequences with the annotation database to get the identification results at different taxonomic levels, including kingdom, phylum, class, order, family, genus, and species. The results are shown in Table 6.

Table 6. Species abundance analysis results

| Sample ID | OTU number | OTU tag number | Kingdom | Phylum | Class | Order | Family | Genus | Species |
|-----------|------------|----------------|---------|--------|-------|-------|--------|-------|---------|
| D3-1 | 37 | 41538 | 1 | 1 | 1 | 3 | 7 | 17 | 26 |
| D3-2 | 43 | 40749 | 1 | 1 | 1 | 4 | 8 | 18 | 27 |
| D3-3 | 43 | 41631 | 1 | 1 | 1 | 3 | 7 | 17 | 25 |
| D4-1 | 75 | 40157 | 1 | 1 | 1 | 3 | 6 | 12 | 16 |
| D4-2 | 95 | 43906 | 1 | 1 | 1 | 5 | 9 | 16 | 19 |
| D4-3 | 98 | 42614 | 1 | 1 | 1 | 4 | 7 | 14 | 18 |
| D5-1 | 60 | 39505 | 1 | 1 | 1 | 4 | 7 | 18 | 25 |
| D5-2 | 59 | 40147 | 1 | 1 | 1 | 4 | 9 | 25 | 33 |
| D5-3 | 60 | 38796 | 1 | 1 | 1 | 5 | 8 | 22 | 28 |
| D6-1 | 40 | 38275 | 1 | 1 | 1 | 5 | 9 | 16 | 20 |
| D6-2 | 48 | 38138 | 1 | 1 | 1 | 5 | 8 | 14 | 18 |
| D6-3 | 49 | 36663 | 1 | 1 | 1 | 4 | 7 | 14 | 18 |
| D7-1 | 40 | 41771 | 1 | 1 | 1 | 5 | 8 | 14 | 19 |
| D7-2 | 46 | 41966 | 1 | 1 | 1 | 4 | 8 | 13 | 15 |
| D7-3 | 46 | 41716 | 1 | 1 | 1 | 5 | 8 | 15 | 19 |
| D8-1 | 49 | 40750 | 1 | 1 | 1 | 4 | 7 | 17 | 24 |
| D8-2 | 47 | 41839 | 1 | 1 | 1 | 4 | 7 | 17 | 24 |
| D8-3 | 50 | 42343 | 1 | 1 | 1 | 5 | 8 | 18 | 25 |

To visually display the composition and abundance of the biological communities at species level in the samples, the system selects the top 18 species in terms of abundance at species level for all samples and generates a stacked bar chart. The chart, shown in Figure 5, illustrates that the species composition and abundance are consistent among the three replicates from the same location. Additionally, a total of 59 fish species were identified across the six locations.




Figure 5. Species composition and abundance stacked histogram

The rank-abundance curve is another function of the system analyzing diversity. It reflects the diversity and evenness of species in the sample by arranging species in descending order of abundance and calculating the relative abundance and ranking of each species. In the horizontal direction, the width of the curve reflects the abundance of species, with a wider range on the x-axis indicating higher species abundance. The shape (smoothness) of the curve reflects the evenness of species distribution in the sample, with a smoother curve indicating a more even distribution of species. The rank-abundance curves for all samples are shown in Figure 6, where the x-axis represents the number of OTUs after ranking, and the y-axis represents the relative abundance of each OTU.




Figure 6. The rank-abundance curve for samples

To display species identification and abundance information at different taxonomic levels (phylum, class, order, family, genus, species) for individual samples more intuitively, the system drawn Krona charts. Taking sample D3 as an example, the Krona chart result is shown below in Figure 7. The circles represent different taxonomic levels from the inner to the outer, and the size of the sectors represents the relative proportion of different identification results.




Figure 7. Krona chart for sample D3

To evaluate whether the sequencing depth can detect most of species, the system drawn Shannon-Wiener curves using the sequencing depth of each sample at different sequencing depths. These curves reflect the biodiversity of each sample at different sequencing depths. When the curve becomes flat, indicates the sequencing data is enough to detect most of the species information in the sample. As shown in Figure 8, we can see the samples curves tend to flatten at 30,000 reads. This suggests that trimming to 50,000 reads is sufficient to detect the majority of species in the samples.




Figure 8. Shannon-Wiener curves of samples

To analyze the differences in species diversity between different locations, the system computed the dissimilarity index of species composition between different samples to assess the degree of species composition differences. The compositional distance analysis plot for all samples is shown in Figure 9, where both the x-axis and y-axis represent the samples, and the color reflects the distance between samples horizontally and vertically. We can see samples from locations D3 and D8 have some similarity in fish species composition, while there is a significant difference in species composition between samples from other locations.




Figure 9. Inter-group distance analysis

Principal Component Analysis (PCA) is another function available in the system for analyzing the differences in species diversity between different locations. The analysis results are shown in Figure 10, where each point represents a sample, and points of the same colour belong to the same group. The closer the distance between two points, the smaller the difference between them. We can see the samples from locations D3 and D8 are closer in distance, indicating a certain similarity in fish species composition between these two locations. This is consistent with the conclusion from the inter-group distance analysis mentioned above.



Figure 10. PCA chart

The OTU sequences were aligned and annotated according to the database to obtain species identification results at different taxonomic levels. The results are shown in Table 8.

Table 8. OTU classification statistics table

| Sample ID | OTU number | OTU tag number | Kingdom | Phylum | Class | Order | Family | Genus | Species |
|-----------|------------|----------------|---------|--------|-------|-------|--------|-------|---------|
| DJ1-1 | 278 | 34884 | 2 | 12 | 20 | 35 | 67 | 102 | 117 |
| DJ1-2 | 316 | 33828 | 2 | 13 | 20 | 34 | 72 | 108 | 123 |
| DJ1-3 | 314 | 33731 | 2 | 13 | 22 | 34 | 75 | 111 | 128 |
| S1-1 | 157 | 36142 | 2 | 9 | 13 | 23 | 42 | 67 | 81 |
| S1-2 | 78 | 35914 | 2 | 8 | 12 | 20 | 32 | 46 | 51 |
| S1-3 | 129 | 34254 | 2 | 10 | 12 | 20 | 39 | 61 | 72 |
| JX1-1 | 1087 | 33282 | 2 | 17 | 27 | 52 | 114 | 206 | 246 |
| JX1-2 | 987 | 36449 | 2 | 14 | 23 | 46 | 110 | 187 | 230 |
| JX1-3 | 1197 | 34670 | 2 | 15 | 25 | 41 | 104 | 195 | 245 |

Sample DJ1, S1 and JX1 are sediment samples collected from Diaojiang, Jiuxu River, and Jiuzhou River in Guangxi. Each sampling location was replicated 3 times. After DNA extraction, COI libraries were prepared using the MGI's solution. The libraries were then subjected to DNBSEQ-G99 sequencing and analysis. The results of the quality control and data processing are shown in Table 7, all samples' Q30 > 96%. Each sample was trimmed to 50,000 reads and conducted to low-quality filtering, denoising, merging, and chimera removal to generate feature sequences.

Table 7. Quality control and data processing results

| Sample ID | Raw reads | Target reads | Down sampled reads | Filtered reads | Q30(%) | Denoised reads | Merged reads | Chimeric reads | Feature reads |
|-----------|-----------|--------------|--------------------|----------------|--------|----------------|--------------|----------------|---------------|
| DJ1-1 | 873230 | 61516 | 50000 | 38417 | 97.34 | 37127 | 35922 | 1038 | 34884 |
| DJ1-2 | 937803 | 61632 | 50000 | 36734 | 96.66 | 35879 | 35112 | 1284 | 33828 |
| DJ1-3 | 824197 | 61491 | 50000 | 36386 | 96.65 | 34886 | 33869 | 138 | 33731 |
| S1-1 | 797266 | 61790 | 50000 | 39331 | 97.56 | 38459 | 37514 | 1372 | 36142 |
| S1-2 | 779515 | 61741 | 50000 | 38541 | 97.24 | 38110 | 37634 | 1720 | 35914 |
| S1-3 | 275573 | 61804 | 50000 | 38650 | 97.57 | 38024 | 36518 | 2264 | 34254 |
| JX1-1 | 253806 | 61694 | 50000 | 38895 | 97.38 | 36448 | 33462 | 180 | 33282 |
| JX1-2 | 352053 | 61795 | 50000 | 40172 | 97.54 | 38325 | 36595 | 145 | 36450 |
| JX1-3 | 205548 | 61814 | 50000 | 38918 | 97.44 | 36917 | 34859 | 189 | 34670 |

Figure 11 shows the stacked bar chart of species composition and abundance at the species level for the three groups of sediment samples. It shows that the species composition and abundance of benthic or planktonic organisms detected in the three replicates of the same sampling location are similar and have similar abundances. However, there are significant differences in species composition and abundance between the three sampling locations.




Figure 11. Species composition and abundance stacked histogram

According to the species identification results of OTUs at different taxonomic levels. Krona plot drawn, as shown in Figure 12. there are significant differences in species composition between the three sampling locations, which is consistent with the above conclusion.



Figure 12. Krona chart for 3 set samples

Three soil samples (T1-T3) collected from the same location were prepared into 18S V4 and ITS1 libraries respectively. The libraries were then sequenced and analyzed on the DNBSEQ-G99. The results of quality control and data processing are shown in Table 9. all samples Q30 > 93%. After trimming to 50,000 reads per sample, the data conducted to low-quality filtering, denoising, merging, and chimera removal. The resulting OTU sequences were then aligned and annotated, as shown in Tables 9 and 10.

In summary, MGI environmental DNA metabarcoding sequencing package has excellent performance, various bioinformatics functions, and compatibility with main types of environmental samples. It is an ideal tool for ecological environmental monitoring and environmental science research.

Table 9. Quality control and data processing results

| Sample ID | Panel | Raw reads | Target reads | Down sampled reads | Filtered reads | Q30 (%) | Denosed reads | Merged reads | Chimeric reads | Feature reads |
|-----------|--------|-----------|--------------|--------------------|----------------|---------|---------------|--------------|----------------|---------------|
| T1 | 18S V4 | 131907 | 64259 | 50000 | 31730 | 96.79 | 30345 | 24598 | 339 | 24259 |
| T2 | 18S V4 | 128451 | 64147 | 50000 | 32647 | 96.91 | 31313 | 24752 | 705 | 24047 |
| T3 | 18S V4 | 93458 | 63723 | 50000 | 33346 | 97.03 | 32071 | 26233 | 404 | 25829 |
| T1 | ITS1 | 188746 | 64919 | 50000 | 47416 | 93.88 | 45077 | 42088 | 507 | 41581 |
| T2 | ITS1 | 227863 | 64846 | 50000 | 47488 | 93.84 | 45245 | 41634 | 165 | 41469 |
| T3 | ITS1 | 218857 | 64773 | 50000 | 47269 | 93.94 | 45132 | 42317 | 503 | 41814 |

Table 10. OTU classification statistics table

| Sample ID | Panel | OTU number | OTU tag number | Kingdom | Phylum | Class | Order | Family | Genus | Species |
|-----------|--------|------------|----------------|---------|--------|-------|-------|--------|-------|---------|
| T1 | 18S V4 | 439 | 24258 | 1 | 37 | 58 | 83 | 90 | 130 | 171 |
| T2 | 18S V4 | 417 | 24046 | 1 | 37 | 53 | 74 | 81 | 115 | 158 |
| T3 | 18S V4 | 452 | 25828 | 1 | 38 | 58 | 83 | 90 | 134 | 175 |
| T1 | ITS1 | 172 | 41581 | 1 | 6 | 16 | 32 | 58 | 73 | 79 |
| T2 | ITS1 | 158 | 41469 | 1 | 7 | 17 | 33 | 59 | 70 | 75 |
| T3 | ITS1 | 164 | 41814 | 1 | 6 | 16 | 31 | 53 | 65 | 72 |

Based on the species identification results of OTUs at different taxonomic levels. Krona plot drawn, as shown in Figure 13.




Figure 13. Krona chart for 2 set soil samples

Ordering information

| Name | Specification | PN |
|---|---------------------------|--------------------|
| Instrument | | |
| DNBSEQ-G99ARS | Include server | 900-000609-00 |
| DNBSEQ-E25RS | Standard Config | 900-000537-00 |
| Automated Nucleic Acid Extractor MGISP-NE32RS | / | 950-000020-00 |
| MGISP-960RS High-throughput Automated Sample Preparation System | Custom Configuration 9-V7 | 900-000154-00 |
| Reagent | | |
| MGIEasy Stool Microbiome Extraction Kit | 48 Preps/Kit | 940-000122-00 |
| MGIEasy Stool Microbiome Extraction Kit | 192 Preps/Kit | 940-000123-00 |
| MGIEasy Tissue Grinding Beads | 96 Preps/ Bottle | 940-000136-00 |
| ATOPlex ITS1 rDNA Library Prep Set | 96 RXN | Group010-000013-00 |
| ATOPlex ITS1 rDNA Library Prep Set | 576 RXN | Group010-000014-00 |
| ATOPlex 18S V4 rDNA Library Prep Set | 96 RXN | Group010-000015-00 |
| ATOPlex 18S V4 rDNA Library Prep Set | 576 RXN | Group010-000016-00 |
| ATOPlex COI mtDNA Library Prep Set | 96 RXN | Group010-000017-00 |
| ATOPlex COI mtDNA Library Prep Set | 576 RXN | Group010-000018-00 |
| ATOPlex Ac12S mtDNA Library Prep Set | 96 RXN | Group010-000019-00 |
| ATOPlex Ac12S mtDNA Library Prep Set | 576 RXN | Group010-000020-00 |
| ATOPlex MiFish Library Prep Set | 96 RXN | Group010-000021-00 |
| ATOPlex MiFish Library Prep Set | 576 RXN | Group010-000022-00 |
| ATOPlex 16SV3V4 rDNA Library Preparation Set | 96 RXN | 940-001261-00 |
| ATOPlex 16SV3V4 rDNA Library Preparation Set | 576 RXN | 940-000725-00 |
| ATOPlex MiFish Panel | 96RXN/Kit | 940-001539-00 |
| ATOPlex MiFish Panel | 576RXN/Kit | 940-001538-00 |
| ATOPlex 16SV3V4 rDNA panel | 96RXN/Kit | 940-001260-00 |
| ATOPlex 16SV3V4 rDNA panel | 576RXN/Kit | 940-000724-00 |
| ATOPlex Ac12S mtDNA Panel | 96RXN/Kit | 940-001535-00 |
| ATOPlex Ac12S mtDNA Panel | 576RXN/Kit | 940-001540-00 |
| ATOPlex COI mtDNA Panel | 96RXN/Kit | 940-001534-00 |
| ATOPlex COI mtDNA Panel | 576RXN/Kit | 940-001537-00 |

| Name | Specification | PN |
|---|------------------|---------------|
| Reagent | | |
| ATOPlex 18SV4 rDNA Panel | 96RXN/Kit | 940-001533-00 |
| ATOPlex 18SV4 rDNA Panel | 576RXN/Kit | 940-001541-00 |
| ATOPlex ITS1 rDNA Panel | 96RXN/Kit | 940-001536-00 |
| ATOPlex ITS1 rDNA Panel | 576RXN/Kit | 940-001532-00 |
| ATOPlex DNA Dual BC Library Prep Set | 96RXN | 940-001191-00 |
| ATOPlex DNA Dual BC Library Prep Set | 576RXN | 940-001190-00 |
| MGEasy Dual Barcode Circularization Kit V1.0 | 16RXN/Kit | 1000020570 |
| DNBSEQ OneStep DNB Make Reagent Kit | 4 RXN | 1000026466 |
| ATOPlex E450 Dual Barcode Balanced Library Reagent | 26ng/Tube | 940-000637-00 |
| Standard Library Reagent (PCR Product) | 1500ng/Tube | 1000027585 |
| High-throughput Sequencing Set | G99 SM FCL PE150 | 940-000410-00 |
| High-throughput Sequencing Set | G99 SM FCL PE300 | 940-000415-00 |
| DNBSEQ-E25RS High-throughput Sequencing Set | FCL PE150 | 940-000567-00 |
| Software | | |
| Metabarcoding Species Identification Software | / | 970-000417-00 |
| Metabarcoding Species Identification Software Package | 96 reports | 970-000456-00 |
| Platform of microorganisms Fast Identification | / | 900-000393-00 |
| Platform of microorganisms Fast Identification and assembly evolution | / | 900-000399-00 |

MGI Tech Co.,Ltd.

📍 Building 11, Beishan Industrial Zone, Yantian District, Shenzhen, China

☎ +86-4000-688-114

🌐 en.mgi-tech.com

✉ MGI-service@mgi-tech.com

